

PACO12727

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGIF2LX

Uniprot

Q8IU1

Synonyms:

TGFB-induced factor homeobox 2-like, X-linked; TGIF2LX; MGC34726; TGIFLX ;

Immunogen:

Human TGIF2LX.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A