

# THOC6 Antibody



PACO12742

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

THOC6

**Uniprot**

Q86W42

**Synonyms:**

THO complex 6 homolog (Drosophila); THOC6; MGC2655; WDR58; fSAP35 ;

**Immunogen:**

Human THOC6.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A