

THRA Antibody



PACO12744

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

THRA

Uniprot

P10827

Synonyms:

thyroid hormone receptor, alpha (erythroblastic leukemia viral (v-erb-a) oncogene homolog, avian); THRA; AR7; EAR7; ERB-T-1; ERBA; ERBA1; MGC000261; MGC43240; NR1A1; THRA1; THRA2; c-ERBA-1 ;

Immunogen:

Human THRA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A