

THTPA Antibody



PACO12745

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

THTPA

Uniprot

Q9BU02

Synonyms:

thiamine triphosphatase; THTPA; MGC2652; THTP; THTPASE ;

Immunogen:

Human THTPA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A