

THYN1 Antibody



PACO12748

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

THYN1

Uniprot

Q9P016

Synonyms:

thymocyte nuclear protein 1; THYN1; HSPC144; MDS012; MGC12187; MY105; THY28; THY28KD ;

Immunogen:

Human THYN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A