

TIGD1 Antibody



PACO12753

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TIGD1

Uniprot

Q96MW7

Synonyms:

tigger transposable element derived 1; TIGD1; EEYORE ;

Immunogen:

Human TIGD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A