

TIGD5 Antibody



PACO12754

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TIGD5

Uniprot

Q53EQ6

Synonyms:

tigger transposable element derived 5; TIGD5; FLJ14926; MGC44883 ;

Immunogen:

Human TIGD5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A