

TIMM23 Antibody



PACO12760

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TIMM23

Uniprot

O14925

Synonyms:

translocase of inner mitochondrial membrane 23 homolog (yeast); TIMM23; MGC22767; PRO1197; TIM23 ;

Immunogen:

Human TIMM23.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A