

TINAGL1 Antibody



PACO12767

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TINAGL1

Uniprot

Q9GZM7

Synonyms:

tubulointerstitial nephritis antigen-like 1; TINAGL1; ARG1; LCN7; LIECG3; TINAGRP ;

Immunogen:

Human TINAGL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A