

# TK1 Antibody



PACO12772

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TK1

**Uniprot**

P04183

**Synonyms:**

thymidine kinase 1, soluble; TK1; TK2 ;

**Immunogen:**

Human TK1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A