

TMCO6 Antibody



PACO12785

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TMCO6

Uniprot

Q96DC7

Synonyms:

transmembrane and coiled-coil domains 6; TMCO6; FLJ39769; PRO1580 ;

Immunogen:

Human TMCO6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A