

TMED10 Antibody



PACO12786

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TMED10

Uniprot

P49755

Synonyms:

transmembrane emp24-like trafficking protein 10 (yeast); TMED10; P24(DELTA); S311125; S3111125; TMP21; Tmp-21-I ;

Immunogen:

Human TMED10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A