

TMED2 Antibody



PACO12787

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TMED2

Uniprot

Q15363

Synonyms:

transmembrane emp24 domain trafficking protein 2; TMED2; FLJ21323; P24A; RNP24 ;

Immunogen:

Human TMED2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A