

TMED4 Antibody



PACO12789

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMED4

Uniprot

Q7Z7H5

Synonyms:

transmembrane emp24 protein transport domain containing 4; TMED4; ERS25; HNFLF ;

Immunogen:

Human TMED4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A