

TMEM106B Antibody



PACO12792

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TMEM106B

Uniprot

Q9NUM4

Synonyms:

transmembrane protein 106B; TMEM106B; FLJ11273; MGC33727 ;

Immunogen:

Human TMEM106B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A