

TMEM108 Antibody



PACO12793

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMEM108

Uniprot

Q6UXF1

Synonyms:

transmembrane protein 108; TMEM108; MGC3040 ;

Immunogen:

Human TMEM108.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A