

TMEM141 Antibody



PACO12796

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TMEM141

Uniprot

Q96I45

Synonyms:

transmembrane protein 141; TMEM141; MGC14141; RP11-216L13.7 ;

Immunogen:

Human TMEM141.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A