

TMEM176A Antibody



PACO12800

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TMEM176A

Uniprot

Q96HP8

Synonyms:

transmembrane protein 176A; TMEM176A; GS188; HCA112 ;

Immunogen:

Human TMEM176A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A