

TMEM41A Antibody



PACO12804

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMEM41A

Uniprot

Q96HV5

Synonyms:

transmembrane protein 41A; TMEM41A;2900010K02Rik; DKFZp686N06149; DKFZp686N107; MGC15397 ;

Immunogen:

Human TMEM41A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A