

TMOD2 Antibody



PACO12810

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TMOD2

Uniprot

Q9NZR1

Synonyms:

tropomodulin 2 (neuronal); TMOD2; MGC39481; N-TMOD; NTMOD ;

Immunogen:

Human TMOD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A