

TMX4 Antibody



PACO12818

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMX4

Uniprot

Q9H1E5

Synonyms:

thioredoxin-related transmembrane protein 4; TMX4; DJ971N18.2; KIAA1162; PDIA14; TXNDC13 ;

Immunogen:

Human TMX4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A