

TNNC1 Antibody



PACO12833

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TNNC1

Uniprot

P63316

Synonyms:

troponin C type 1 (slow); TNNC1; CMD1Z; TNC; TNNC ;

Immunogen:

Human TNNC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A