

TNNI1 Antibody



PACO12834

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TNNI1

Uniprot

P19237

Synonyms:

troponin I type 1 (skeletal, slow); TNNI1; DKFZp451O223; SSTNI; TNN1 ;

Immunogen:

Human TNNI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A