

TNNT2 Antibody



PACO12837

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TNNT2

Uniprot

P45379

Synonyms:

troponin T type 2 (cardiac); TNNT2; CMH2; CMPD2; MGC3889; RCM3; TnTC; cTnT ;

Immunogen:

Human TNNT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A