

# TNS3 Antibody



PACO12840

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TNS3

**Uniprot**

Q68CZ2

**Synonyms:**

tensin 3; TNS3; DKFZp686K12123; DKFZp686M1045; FLJ13732; FLJ35545; H\_NH049123.2; MGC88434; TEM6; TENS1 ;

**Immunogen:**

Human TNS3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A