

# TOB1 Antibody



PACO12842

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TOB1

**Uniprot**

P50616

**Synonyms:**

transducer of ERBB2, 1; TOB1; APRO6; MGC104792; MGC34446; PIG49; TOB; TROB; TROB1 ;

**Immunogen:**

Human TOB1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A