

TOLLIP Antibody



PACO12844

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TOLLIP

Uniprot

Q9H0E2

Synonyms:

toll interacting protein; TOLLIP; FLJ33531; IL-1RAcPIP ;

Immunogen:

Human TOLLIP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A