

TOMM40 Antibody



PACO12849

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TOMM40

Uniprot

O96008

Synonyms:

translocase of outer mitochondrial membrane 40 homolog (yeast); TOMM40; C19orf1; D19S1177E; PER-EC1; PEREC1; TOM40 ;

Immunogen:

Human TOMM40.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A