

# TOR2A Antibody



PACO12855

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TOR2A

**Uniprot**

Q8N2E6

**Synonyms:**

torsin family 2, member A; TOR2A; FLJ14771; MGC99558; TORP1 ;

**Immunogen:**

Human TOR2A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A