

TPPP2 Antibody



PACO12871

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TPPP2

Uniprot

P59282

Synonyms:

tubulin polymerization-promoting protein family member 2; TPPP2; C14orf8; P18; p25beta ; TPPP/p18; Protein p25-beta

Immunogen:

Human TPPP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A