

TRADD Antibody



PACO12881

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRADD

Uniprot

Q15628

Synonyms:

TNFRSF1A-associated via death domain; TRADD; Hs.89862; MGC11078 ;

Immunogen:

Human TRADD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A