

TRAF7 Antibody



PACO12887

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRAF7

Uniprot

Q6Q0C0

Synonyms:

TNF receptor-associated factor 7; TRAF7; DKFZp586I021; MGC7807; RFWD1; RNF119 ;

Immunogen:

Human TRAF7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A