

TRAPPC2 Antibody



PACO12893

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRAPPC2

Uniprot

P0DI81

Synonyms:

trafficking protein particle complex 2; TRAPPC2; MIP-2A; SEDL; SEDT; TRS20; ZNF547L; hYP38334 ;

Immunogen:

Human TRAPPC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A