

TRAPPC4 Antibody



PACO12895

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRAPPC4

Uniprot

Q9Y296

Synonyms:

trafficking protein particle complex 4; TRAPPC4; CGI-104; HSPC172; PTD009; SBDN; SYNBINDIN; TRS23 ;

Immunogen:

Human TRAPPC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A