

TREM1 Antibody



PACO12898

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TREM1

Uniprot

Q9NP99

Synonyms:

triggering receptor expressed on myeloid cells 1; TREM1; TREM-1 ;

Immunogen:

Human TREM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A