

# TRERF1 Antibody



PACO12899

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRERF1

**Uniprot**

Q96PN7

**Synonyms:**

transcriptional regulating factor 1; TRERF1; BCAR2; FLJ21198; HSA277276; RAPA; RP1-139D8.5; TReP-132; dJ139D8.5 ;

**Immunogen:**

Human TRERF1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A