

TRIM24 Antibody



PACO12909

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRIM24

Uniprot

O15164

Synonyms:

tripartite motif-containing 24; TRIM24; PTC6; RNF82; TF1A; TIF1; TIF1A; TIF1ALPHA; hTIF1 ;

Immunogen:

Human TRIM24.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A