

TRIM46 Antibody



PACO12920

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TRIM46

Uniprot

Q7Z4K8

Synonyms:

tripartite motif-containing 46; TRIM46; FLJ23229; GENEY; TRIFIC ;

Immunogen:

Human TRIM46.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A