

TRIM6 Antibody



PACO12923

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRIM6

Uniprot

Q9C030

Synonyms:

tripartite motif-containing 6; TRIM6; RNF89 ;

Immunogen:

Human TRIM6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A