

TRMT1 Antibody



PACO12930

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRMT1

Uniprot

Q9NXH9

Synonyms:

TRM1 tRNA methyltransferase 1 homolog (S. cerevisiae); TRMT1; FLJ20244; TRM1 ;

Immunogen:

Human TRMT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A