

TRMT5 Antibody



PACO12933

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRMT5

Uniprot

Q32P41

Synonyms:

TRM5 tRNA methyltransferase 5 homolog (S. cerevisiae); TRMT5; KIAA1393; MGC111453; TRM5 ;

Immunogen:

Human TRMT5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A