

TRPC6 Antibody



PACO12940

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRPC6

Uniprot

Q9Y210

Synonyms:

transient receptor potential cation channel, subfamily C, member 6; TRPC6; FLJ11098; FLJ14863; FSGS2; TRP6 ;

Immunogen:

Human TRPC6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A