

# TRPV4 Antibody



PACO12943

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRPV4

**Uniprot**

Q9HBA0

**Synonyms:**

transient receptor potential cation channel, subfamily V, member 4; TRPV4; OTRPC4; TRP12; VR-OAC; VRL-2; VRL2; VROAC ;

**Immunogen:**

Human TRPV4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A