

TSC22D3 Antibody



PACO12948

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSC22D3

Uniprot

Q99576

Synonyms:

TSC22 domain family, member 3; TSC22D3; DIP; DKFZp313A1123; DSIPI; GILZ; TSC-22R; hDIP ;

Immunogen:

Human TSC22D3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A