

# TSC22D4 Antibody



PACO12949

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

TSC22D4

**Uniprot**

Q9Y3Q8

**Synonyms:**

TSC22 domain family, member 4; TSC22D4; THG-1; THG1 ;

**Immunogen:**

Human TSC22D4.

**Storage:**

0.1M NaHCO<sub>3</sub>, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A