

# TSPAN15 Antibody



PACO12960

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSPAN15

**Uniprot**

O95858

**Synonyms:**

tetraspanin 15; TSPAN15;2700063A19Rik; NET-7; TM4SF15 ;

**Immunogen:**

Human TSPAN15.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A