

TSPAN6 Antibody



PACO12963

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSPAN6

Uniprot

O43657

Synonyms:

tetraspanin 6; TSPAN6; T245; TM4SF6; TSPAN-6 ;

Immunogen:

Human TSPAN6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A