

# TSR1 Antibody



PACO12966

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSR1

**Uniprot**

Q2NL82

**Synonyms:**

TSR1, 20S rRNA accumulation, homolog (S. cerevisiae); TSR1; FLJ10534; KIAA1401; MGC131829 ;

**Immunogen:**

Human TSR1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A