

TTC39A Antibody



PACO12972

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TTC39A

Uniprot

Q5SRH9

Synonyms:

tetratricopeptide repeat domain 39A; TTC39A; C1orf34; DEME-6; KIAA0452 ;

Immunogen:

Human TTC39A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A