

TTC4 Antibody



PACO12973

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TTC4

Uniprot

O95801

Synonyms:

tetratricopeptide repeat domain 4; TTC4; DKFZp781B0622; FLJ41930; MGC5097 ;

Immunogen:

Human TTC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A